SAM[®] NEO

MANUAL

Anatomically Correct Simulation

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
The **Upper Member** of the articulator consists of **Condylar Housings** posteriorly and **Upper Incisal Guidance Table** anteriorly. Together, they simulate the TMJ and upper anterior teeth guidance.

The moving **Lower member** of the articulator with **Condylar Elements** posteriorly and the **Incisal Pin** anteriorly, together, they simulate the mandibular condylar movement and lower anterior teeth guidance.

This, therefore, provides for a very logical and natural laboratory simulation of the patient movements which are very easy to understand and reproduce.

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The red dot in the center of the lower member condylar element corresponds to a hinging axis rotation center of the mandibular condyles. Also known as "Hinge Axis".

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Above the red dotted/dashed line is the axis-orbital plane (AOP)

The upper anterior incisal table surface along with the posterior condylar element centers of rotation are on the same horizontal plane.

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
The lower incisal pin in contact with the upper incisal guidance table is used to simulate the relationship of upper and lower incisal edges of the patient's teeth. Thus promoting normal functional patterns.

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
The SAM® Articulator System


All current articulators and instruments made by SAM® are based on the height „P“ (prosthetic height = 121 mm between upper and lower members). The intercondylar distance is 110 mm. A particular feature of the system is the SAM® 15 mm distal offset arrangement of the mounting plates in the articulator. This creates a more normal patient simulation of mounted models in the SAM® articulator.

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The SAM® NEO and all articulators made by SAM® are ARCON Instruments. To simulate the patient's authentic TMJ movements, the lower articulator member condylar ball moves in the upper member condylar housing and this movement corresponds to the TMJ anatomical model. This natural anatomical movement is logically simulated and understandable.

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


Caution: Do not touch the end of the spring bolt (as shown here) with your thumb!

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The centric lock is open when the spring is compressed. The position of the bolt can be fixed with the centric lock screw.

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When the locking screw is open, the pressure spring simulates the muscular strength of the patient's jaw. The spring resistance repeatedly moves the condyle ball into the centric position.

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On the other hand, if the centric lock is completely open, sliding movements can be performed.

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Both Centric Locking Assemblies must be OPEN in order to make Articulator Movements.

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often called as **Bennett angle** is a Mediotrusion curve and no straight pathway, as in most mean value articulators. The Bennett guide is already given in the condylar housing and corresponds to the most common patient situation in a healthy TMJ.

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





pathway
顆路

The Bennett movement is a downward, inward, forward movement or downwards, forward, inward movement all at the same time of the condyle on the mediotrusion (bal-ancing, non-working) side while the laterotrusion (working, rotating) side moves virtually straight lateral for about 0.3 mm.

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An exaggerated mediotrusion movement is shown here. The upper and lower jaws move in the opposite direction. The canine guidance causes disclusion of the posterior teeth when the mandible moves laterally and/or in a pro-trusive direction.

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To perform a laterotrusion movement, the locking screw of the spring-loaded bolt is loosened.

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By holding the condylar element firmly against the medial wall of the condylar housing and pushing the condylar element downward, inward, and forward the condylar element guiding on the condylar housing pathway will move the lower member towards the midline to perform the mediotrusion movement.

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Condylar Pathway Inclination

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In the SAM® NEO Articulator you can adjust the condylar inclination pathway. The condylar pathway curvature is the equivalent of the SAM® white curva-ture pathway.

The joint path is most often a curved trajectory. This was reflected in the condylar SAM® NEO.

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
Setting the Condylar Pathway Inclination. The adjustment of the inclination is performed by loosening the posterior Thumb screw on the rear of the upper member of the articulator.

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
The inclination of the Condylar Pathway is adjustable from 15° to 65°.

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The path shown is the real path of the center of rotation of the temporomandibular joint head of the patient.

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To simulate this as faithfully as possible, the condyle of the envelope of condylar pathway elements (joint path curvature) must follow.

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
Insert the Tilt Support Rods into the holes on the back of the articulator posts. They support the articulator in tilted positions in an angle of approx. 40°.

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This improves the view of the anterior teeth. Occasionally apply silicone spray, to maintain the O ring (ART 526002).

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The SAM® NEO is like all SAM® articulators which are designed originally for an upside down position as well as upright. Notice how this improves the look at the occlusion while observing and controlling the function and effect of the condylar elements of the articulator.

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Insert the vertical support rod into the hole behind the thumb screw of the incisal table.

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It supports the upper member of the articulator when the articulator is open.

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This improves the view of the occlusion, and makes it easier to work on the model (e.g wax up).

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Occasionally apply silicone spray to maintain the support rod and the inner O ring (ART 525002).

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The SAM articulators are available with

1. Screw Type Mounting Plate System (ART120)
2. MPS Magnetic Plate System (MPS150)
3. AXIOSPLIT® System Plate Oval (ASP420)
4. AXIOSPLIT® system Plate Round (ASP425)

The different SAM® Mounting Plate Systems are available for all SAM® articulators and equipment. All SAM® instruments can be converted without difficulty.

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Mounting Plate with Screw Type System

The yellow mounting plate is secured with a Thumb Screw and this provides maximum retention for assembly of models with small amounts of plaster. The mounting plates are intended for single use and for storage and easy reattachment of the mounted models.

The low cost of the mounting plates indicates this is the procedure of choice. Not to mention the saving of the costly time to try and recycle the plates. The positioning of dental cast in the articulator is carried out very precisely by means of two precision alignment pins.

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Benefit: the patient model remains from beginning to end in the same articulator

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The SAM® Magnetic Mounting Plate MPS

The MPS Magnetic Mounting Plate System is an alternative to the historical screw mounting plate.

This Magnetic Plate System allows for the quick removal and replacement of the models in the articulator. The conversion of screw type to magnetic system can be done easily and quickly with the appropriate kits (MSP 101K and 101K MPS).

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Benefit: Quick exchange of the casts

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The AXIOSPLIT® system

The AXIOSPLIT® system is used for trouble-free and accurate exchange of models between different articulators (for example, old and new devices), in practice and laboratory. So the models need only to be sent. The articulators stay in place.

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A Master Check ensures the synchronizing with all SAM® devices. The mounted models can be interchanged without loss of fitting accuracy in all SAM® devices. You will initially need the AXIOSPLIT® Starter Kit with Master Check, norm incisal pin and related matrix holder and an articulator kit for each additional instrument. The exact procedure and more information is described in our AXIOSPLIT-Manual.

Benefit: Highest precision and your own synchronizing

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The Screw Type Mounting Plate has a round hole behind the threaded nut and an elongated slot anteriorly to assure precise plate positioning on all SAM® instrumentation. set to the precision mounting plate pins. The elongated hole shows to the front.

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
With one hand the plate is held in position and with the second hand thumb and fore finger the thumb screw is tightened to contact. It is further tightened firmly only 1/4 of a turn. tightened with the knurled screw on the other side.

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The mounting plate is fully seated when it sets firmly on the border areas without gap.

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Take note: The mounting plate is screwed firmly to the articulator, the mounted model may not be set easily and quickly removed on and off as like with the MPS system.

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Note: The mounting plate is disposable and should not be used more than once since the difficulty of removing set gypsum and the cleaning of the mounting plate normally will destroy plate and destroy its accuracy.

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MPS-Magnetic Mounting Plates Attachment MPSタイプ

The adhesion plate is set with the rounded edge downwards into the provided hole.

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
The MPS-mounting plate is attached with the round side containing the slot anteriorly onto the positioning pins.

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
Take note: The slightly rounded and shiny side shows downwards when it is inserted into the mounting plate (the rounded edge side shows space for plaster).

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To reproduce or build tooth guidances there should to be an incisal table. It must be able to produce the inclination and direction of the tooth guides.

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
To simulate the anterior teeth and canine guidance authentically, the incisal table is always attached to the upper part of the articulator.

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The Plane Incisal Table

The incisal table is located as a guiding surface on the top member of the articulator. The flat plane incisal table is used to secure the vertical dimension and to support the incisal pin and can be individualized with plastic via the lateral wings.

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Adjusting the Adjustable Incisal Table(Accessories)

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1 The incisal table is not pushed completely into the articulator top until it stops until the anterior moveable section has been completely raised.

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3 Set incisal table on steepest inclination 80° and the rotation to 0°

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4 The incisal table is seated fully in the slot on the upper member of the articulator until the tip of the incisal pin simultaneously touches the inclined and the horizontal surface of the incisal table. Then the knurled screw is tightened.

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5 In order to check the correct position of the incisal table, the incisal table can be rotated anywhere from 0° to 60° around its own axis when the articulator is closed.

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The incisal pin should be always located in the lower part of the SAM® articulator.

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
The zero mark on the scale in the middle of the incisal pin is normally the starting position for the dental laboratory. The incisal pin can be fixed in the restorations bite height. The incisal pin can be raised to 11 mm or lowered up to 6 mm.

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Due to the articulator construction the desired value must be set twice at the incisal pin while raising the occlusion, so e.g. in an increase of 2.5 mm in the occlusion the adjustment at the incisal pin is +5 mm.

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When lowering of the bite to 1 mm, the incisal pin should be lowered 2 mm.

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The micro adjustable incisal pin can be adjusted accurately and within 0.05 mm range. When lateral knurled screw is opened the pin can be finely adjusted on the upper black scale of incisal continuously.

This can be finely adjusted in the micrometer range.

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Example of raising by 2.5 mm

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Open the lateral knurled screw slightly (about 1/8 turn).

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
Rotate the black micro adjustment screw clockwise until the horizontal mark on the scale is in the center of the pin and shows "2" with the vertical mark on the upper black scale "5".

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The position is fixed with the lateral knurled screw.

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AXIOQUICK®

The axis orbital plane in the patient corresponds in the articulator to the line between the bearing surface of the incisal table to the incisal pin (1) and the centers of the condylar balls (hinge axis) (2).

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
*Cranium-related mounted models in the articulator are an essential part of diagnosis and therapy.
The correct position of the maxillary model to the reference plane is an absolute prerequisite for model mounting with a centric registration, but also for setting condylar pathways and even when using mean values.*

*The Frankfort Horizontal (FH) Plane represents a reference plane for the patient and the articulator.
The scissors mechanism of the SAM® anatomic facebow creates and allow for an optimal positioning of the transfer fork (from a functional and aesthetic point of view).*

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*Please find the exact procedure and the correct technique to attach the facebow and the proper maintenance and care in our AXIOQUICK® facebow manual.
This manual can be provided on request by SAM® or can be downloaded as PDF file.*

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www.sam-dental.de/downloads


Model mounting in the articulator

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The transfer stand SE is used for an easy and fast maxilla mounting in all SAM® articulators and laboratory cast mounting instruments. After the transfer the anatomic facebow can remain in the treatment room, only the transferfork assembly is needed for the model transfer. All versions of the SAM® transfer fork assemblies are compatible with the transfer stand AX. With the transfer stand the upper jaw position of the patient can be quickly and precisely transferred to the articulator.

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A more detailed manual for model mounting can be provided on request by SAM® or can be downloaded as PDF file.

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The occlusal plane indicator set (ART 385) with a plane and a curved plate with Speescher curve can be used for the average-maxillary model mounting or for a prosthetic set up.

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Flat occlusal plane indicator 平面タイプ

The height of the stage device is set to 55 mm (corresponding to the halving of the height - 110 mm between upper and lower part). The medium order setting the depth is 75 mm.

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Average Upper Jaw Model Mounting

The height of the stage device is set to 55 mm (corresponding to the halving of the height - 110 mm between upper and lower part) The medium order setting the depth is 75 mm .

The angle scale can be adjusted between 0 ° - 45 °, we recommend that the inclination is usually between 0° - 15°.

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
Curved Template 曲面タイプ


Curved Template is used in the production of full dentures and as an aid for a setting up of the maxillary teeth. It provides a Spee Curve. The marks on the template indicate auxiliary lines as the middle of jaw ridge. In the loosened state, the cap may also adjust the inclination angle of the remaining teeth can be determined as occlusion plane.

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Cleaning and Maintenance of the Articulator

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The small centric screw on condylar housing, and the set screws on the incisal pin and table should be occasionally cleaned with mild grease solvents and then lubricated with e.g. Multi-use lubricant (ART 193).

The SAM® NEO is an easy to care for articulator. Plaster residues can be removed with water, soap and a soft brush.

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To remove greasy stains, we recommend mild grease solvents. Condylar ball and housing should always be clean and dry (oil and grease free).

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Accessories and Consumable supplies for SAM® NEO

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ART 114K

Mounting Plate Thumb Screw Assembly

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MPS 101K

MPS additional articulator kit yellow

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ART 533

incisal pin III PM, micro adjustable, complete

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magnetic block for transfer

fork support

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All of our products are all manufactured in Germany.

Made in Germany